



AIR PRO IPL

INTENSE PULSED LIGHT

HAIR REMOVAL DEVICE
USER GUIDE


Thank you for purchasing the GLO AIR PRO IPL!

Please read the manual carefully before use and store the device properly.


The pictures in this manual are for illustration purposes only.

Device Description


1.1. Instrument Description



1.2. IPS Display Interface




1.3. Accessories




How (IPL) Intense Pulsed Light Device Works


The GLO AIR PRO IPL (Intense Pulsed Light) device emits light pulses with a wavelength range of 550nm to 1200nm onto the skin. Melanin in the hair absorbs these pulses, which heats and destroys the hair follicle. This painless process encourages the hairs to enter the resting phase, leading to long-lasting hair-free skin.




Darker hair absorbs the light pulses more effectively.



The light pulses stimulate the hair follicle to enter a resting phase.



The GLO AIR PRO IPL uses intense pulsed light to stimulate the regeneration of collagen and elastin fibers, improving the skin's color and texture without the need for surgery.




The GLO AIR PRO IPL can address visible damage caused by sun exposure, known as photoaging, which often appears on the face, neck, hands, or chest. After using the device, you can expect a more even skin tone, reduced wrinkles, and overall skin rejuvenation.


Contraindications

1. General Conditions

- ❖ Do not use the device on dark skin. Refer to the skin tone table below to determine if your skin tone is suitable for the device.



- ❖ The device is not effective on hair colors such as light brown, light gold, red, or white gray. Refer to the hair color table below to check if your hair color is suitable for using the device.



- ❖ Never use the device if you are pregnant or breast feeding.
- ❖ Never use the device if you have active implants.

Contraindications

2. Medications / History

- ❖ Never use the device if you are taking any of the medications/treatments listed below:
 - I) If your skin is currently being treated with or has been treated in the past week with Alpha-Hydroxy Acids (AHAs), Beta-Hydroxy Acids (BHAs), topical isotretinoin, or azelaic acid.
 - II) If you have taken isotretinoin (Accutane or Roaccutane) in the last 6 months, as this treatment can make your skin more susceptible to tears, wounds, and irritation.
 - III) If you are taking photo-sensitizing agents or medications, check the package insert for any warnings about photo-allergic or photo-toxic reactions, or if you need to avoid sunlight while using the medication.
 - IV) If you take anticoagulants, including heavy use of aspirin, without allowing for a minimum 1-week washout period before each treatment.
 - V) If you have received radiation therapy or chemotherapy within the past 3 months.
 - VI) If you are on painkillers that reduce skin sensitivity to heat.
 - VII) If you take immunosuppressive medications.
 - VIII) If you have or had surgery in the areas to be treated for IPL in the last 3 weeks.

3. Disorders

- ❖ Never use the device;
 - I) If you have diabetes or other systemic or metabolic diseases.
 - II) If you have congestive heart disease.
 - III) If you have a disease related to photo-sensitivity.
 - IV) If you have a history of collagen disorder.
 - V) If your skin is sensitive to light and easily develops a rash or an allergic reaction.
 - VI) If you have a skin disease such as active skin cancer or you have a history of skin cancer or any other localized cancer in the areas to be treated by IPL.
 - VII) If you have a history of vascular disorder.
 - VIII) If you have any bleeding disorder.
 - IX) If you have a history of immunosuppressive disease.

Contraindications

4. Skin Conditions

- ❖ Never use the device;
 - I) If you have infections, eczema, burns, inflammation of hair follicles, open lacerations, herpes simplex (cold sores), wounds or lesions and hematomas in the area to be treated by the IPL.
 - II) Or irritated (red or cut), sunburned, recently tanned or fake-tanned skin.
 - III) On these areas (without consulting your doctor first) such as moles, freckles, large veins, darker pigmented areas, scars and skin anomalies. This can result in a burn and a change in skin color, which makes it potentially harder to identify skin-related diseases.
 - IV) On areas such as; warts, tattoos or permanent makeup.

5. Areas

- ❖ Never use the device;
 - I) Around the eyes or near the eye brows.
 - II) On the lips, nipples, areola, labia minora, vagina, anus and inside of the nostrils and ears.
 - III) Men must not use the device on the face and neck including all beard growing areas, nor on the whole genital area.
 - IV) On areas where you use long lasting deodorants.
 - V) Over or near anything artificial like silicone implants, subcutaneous injection ports (for instance, an insulin dispenser). or piercing.

⚠ If you are not sure whether you can use the device, consult your doctor.

Important Notes

1. Warning

- ❖ Keep the device and the adapter dry.
- ❖ If the device is broken, do not touch any inner part to avoid electric shock.

Important Notes

- ❖ Water and electricity are a dangerous combination. Do not use this device in wet surroundings.
- ❖ This device is not intended for use by individuals (including children) with reduced physical, sensory, or mental capabilities unless they are supervised or instructed on its use by a responsible person.
- ❖ This device is not intended for use by children under the age of 15.
- ❖ Always inspect the device before use. Do not use the device or adapter if it is damaged.
- ❖ Do not use the device, if the flash (light) emission window is broken.
- ❖ The light emission window can become very hot after use. Avoid touching the inner part of the light emission window until it has cooled down.

2. Device Care

- ❖ Ensure that the airflow through the device's ventilation slots is not obstructed.
- ❖ Keep the device stored in a clean, dry location free from dust.
- ❖ To avoid damage, keep the device away from direct sunlight and UV exposure for extended periods.

3. Caution

- ❖ For hygiene purposes, the device should be used by only one person.
- ❖ Use the device only with settings appropriate for your skin type. Using higher settings than recommended may increase the risk of skin reactions and side effects.
- ❖ Do not look directly at the light emission window while using the device. Although protective eye-wear is not mandatory, it is provided with the device for your safety as a precaution.
- ❖ Use the device in a well-lit room to reduce glare and make the light less harsh on your eyes.
- ❖ Never leave the device unattended while it is switched ON. Always turn OFF the device after use.
- ❖ Tanning with natural or artificial sunlight can affect the sensitivity and color of your skin.
- ❖ Before using the device, ensure your skin is clean, hair-free, completely dry, and free of oily substances.
- ❖ Do not treat the same skin area more than once in a single session.
- ❖ The device should never cause pain. If you experience discomfort, lower the light intensity level.

Course of Treatment

1st Treatment - (Week - 0)


- ❖ After the first treatment, it may take 1 to 2 weeks for the hair to fall out.
- ❖ During the initial treatments, you might notice some hair growth. This is likely due to hairs that were in their growth phase during the first treatment.

2-3 Treatments - (Week 2 - Week 4)

- ❖ After 2-3 treatments, you should observe a noticeable reduction in hair growth. To effectively treat all hair, continue following the recommended treatment schedule.

4-5 Treatments - (Week 6 - Week 8)




- ❖ After 4-5 treatments, you should see a significant reduction in hair growth in the areas you treated with the device.
- ❖ A reduction in hair density should also be noticeable. Continue with regular touch-ups every 4 weeks to maintain the results.



Before Using

1. Prepare Your Skin

- ❖ Before using the GLO AIR PRO IPL, remove all visible hair from the area you wish to treat.
- ❖ Ensure there are no hairs remaining above the surface of the skin. This allows the light to be easily absorbed by the melanin in the hair roots, ensuring effective treatment.
- ❖ You should only shave before using the IPL device.
- ❖ Do not wax, epilate, or use depilatory creams, as these methods remove hair from the roots (follicles) and hinder the IPL process.
- ❖ If you choose to wax, wait 6-8 days before using the device to allow your skin to rest.
- ❖ We recommend taking a shower before treatment to ensure your skin is free from oil and dirt, and then pat your skin dry.



2. Skin Test

- ❖ Check the skin tone table to ensure your skin type is suitable for using the device.


3. Hair Test

- ❖ Check the hair color table to verify if your hair color is suitable for using the device.




While Using

1. Selecting the Right Intensity


- ❖ The device offers 8 different light intensities and recommends the appropriate setting based on your skin tone. You can always choose the light intensity that you find most comfortable.



- ❖ Press  for 2 seconds to turn the device ON.
- ❖ Press the  button briefly to adjust the light intensity level. Each  press will increase the intensity gradually from the 1st mode to the 8th mode. After reaching the 8th mode, pressing the  button will decrease the intensity gradually back to the 1st mode.



- ❖ Press  and hold the button for 2 seconds to switch between "AUTO" and "MANUAL" mode.
- ❖ Press  and hold the  button for 2 seconds to switch between "HAIR REMOVAL (HR)" and "SKIN REJUVENATION (SR)" modes.



Note: There is only one intensity level for skin rejuvenation.


While Using

For hair removal, select the desired intensity setting. Refer to the table below to choose the appropriate light intensity setting.


Skin Type & Tone	Light Intensity Setting
I White: always sunburns, never tans.	6-8
II Beige: easily sunburns, tans minimally.	5-7
III Light brown: easily sunburns, tans slowly to a light brown.	4-6
IV Middle brown: rarely sunburns, tans easily.	2-4
V Dark brown: rarely sunburns, tans very easily.	You're unable to use the device.
VI Brownish black or darker: rarely or never sunburns, tans very dark.	You're unable to use the device.

- ❖ The GLO AIR PRO IPL should never cause pain. If you feel uncomfortable, reduce the light intensity setting.
- ❖ After moving to another body area or recent tanning, perform a skin test to determine the appropriate light intensity.


2. Handling the Device



- ❶ Clean the light emission window before use.




- ❷ Insert the plug.




- ❸ Plug the power supply into the wall socket.



- ❹ Switch on the device and select the appropriate light intensity for your skin tone.



- ❺ Place the device flat against your skin at a 90° angle.



- ❻ Ensure the device is flush against your skin with no gaps.


While Using

- ❖ Press the device firmly against your skin to ensure proper contact. Press the  button to release a flash. You may feel a warm sensation, which is normal.
- ❖ Place the device's light window in full contact with the skin to ensure it works correctly. If it's not in proper contact, the device will not function.
- ❖ Move the device to the next target area. After each flash, wait 1-3 seconds until the device is ready for the next flash.
- ❖ For larger body areas, use AUTO mode. Hold down the button and move the device along your skin.
- ❖ For smaller areas, switch to MANUAL mode.
- ❖ When you have finished the treatment, press and hold  for 2 seconds to switch OFF the device.
- ❖ Remove the adapter from the electric socket.
- ❖ It is important to wipe the light window clean after each use.

3. Selecting the Right Intensity

a. Hair Removal Schedule

- ① In the first phase of the treatment, use the device once a week for 4-8 weeks. (Depending on your skin's response).
- ② After this period, use the device once a month for 3 months or until you are satisfied with the results.
- ③ To maintain smooth, hair-free skin, use the device every 2 to 3 months or as needed.



b. Skin Rejuvenation Schedule

For skin rejuvenation, we recommend using the device every 3 days.

After Use

1. After Care

- ❖ After the IPL treatment, you can safely apply lotion, creams, or moisturizer to the treated area.
- ❖ If you experience skin irritation or redness after the treatment, wait until it subsides before applying any products, such as lotion or cream, to your skin.
- ❖ If you experience skin irritation or redness after applying any product, such as lotion or cream, wash it off with water.

2. Common Reactions

- ❖ If your skin shows redness or feels prickly, tingly, or warm, this reaction is harmless and will disappear soon.
- ❖ Dry skin and itching may occur due to shaving or a combination of shaving and IPL treatment. To soothe the area, use an ice pack or a damp, cool cloth. If dryness persists, apply a non-scented moisturizer to the treated area.

3. Rare Side Effects

- ❖ Burns or excessive redness (e.g. around hair follicles).
- ❖ Skin discoloration occurs very rarely.
- ❖ Epidermal heating.

Device Cleaning and Storage

- ❖ After use, switch OFF the device, unplug it, and allow it to cool down.
- ❖ Use a soft, moist cloth to clean the GLO AIR PRO IPL, including the light window, the outside surface of the attachment, and the reflectors on the side.
- ❖ Store the device in a dust-free, dry place at room temperature between -25°C and 75°C.

Device Specifications

Product	IPL Hair Removal Device	Spectrum	>550nm
Model	GLO AIR PRO IPL	Pulse Duration	0.9+/-0.2ms
Rated Input	AC100-240V/50-60Hz	Max Light Energy	7.92-15.31J
Rated Output	DC 12V, 2.5A	Operating Temperature	5-35°C
Executive Standard	GH4706.1 GB4706.15	Operating Humidity	15%-75%
Protection Against Electric Shock	II	Weight	209g
Protections Rating	IP 30(EN 60529)	Size	H179*W70*D35mm

CERTIFICATE

Date of Manufacture _____

Inspection Result _____

qualified
QC002

Inspector No _____



EMC



FC

